MEDIFLOW ULTRA II

LOW PRESSURE REGULATORS

- GB
- HU
- CZ
- PL
- RU

INSTRUCTION FOR USE

Document No.:	735 800 000 096
Date of issue:	01/03/2007
Revision No.:	1

1. Foreword

This instruction for use covers essential information for the whole life cycle of the regulators:

- Operation
- Maintenance
- Cleaning
- Disposal

2. Intended Use

MediFlow Ultra is a low pressure regulator with flow-selector device intended for the administration of the following medical gases in the treatment, management, diagnostic evaluation and care of patients:

- oxygen;
- oxygen-enriched air;
- nitrous oxide;
- medical air;
- synthetic air;
- carbon dioxide;
- helium;
- xenon;
- specified mixtures of the gases listed above

The product is not intended to be use with air or nitrogen for driving surgical tool.

This low pressure regulator is intended to be fitted to medical gas pipeline system terminal units in hospital or ambulance car or to quick connector outlet of medical regulator, with nominal pressures up to 800 kPa.

3. Operational, transport and storage safety requirements

	Operations	Transport	Storage
 <ul style="list-style-type: none"> • Keep away the product and its associated equipments from <ul style="list-style-type: none"> - heat sources (fire, cigarettes, ...) - flammable materials - oil or grease - water - dust 	✓	✓	✓
 <ul style="list-style-type: none"> • The product and its associated equipments must be prevented from falling over. 	✓	✓	✓
 <ul style="list-style-type: none"> • Use only the product and its associated equipments in well vented area. 	✓		
 <ul style="list-style-type: none"> • Always maintain oxygen cleanliness standards. 	✓		

Before initial use the product shall be kept in its original packaging. GCE recommends use of the original packaging (including internal sealing bag and caps) if the product is withdraw from operation (for transport and storage).

Statutory laws, rules and regulations for medical gases, accident prevention and environmental protection must be observed.

4. Personnel instructions and training

According to Medical Devices Directive 93/42/EEC the product owner must ensure that all personnel handling the product are provided with the operating instructions & performance data and are fully trained to carry out that operation. Trainees need to be supervised by an experienced person.

5. Product description

Gas from the medical gas pipeline system entering through the inlet connector, then calibrated in the flow control head and deliver to the user through the flow outlet user connection. Different flows can be selected by the flow control knob on the flow head.

A. Inlet quick connector (QC)

The low pressure regulator is fitted to the medical gas pipeline system terminal unit or to medical regulator by a gas specific male probe.

B. Flow control head

Different flows can be selected by the flow control knob on the flow head. Depending on product variant, the flow control head is fitted with a stop function in between maximum flow and zero position.

C. User outlet

The "flow outlet" user connection is the outlet from the flow control head normally with a push on hose fitting (hose nipple) or threaded (for humidifier).

This outlet is for supplying a gas flow (l/min) at atmospheric pressure directly to a patient, e.g. through a cannula or a facemask..

- **The product is not fitted with a pressure relief valve. The user must ensure that an over pressure protection is in place on the medical gas supply pipeline (see requirement of ISO 7396-1) or on medical regulator (see requirement of ISO 10524-1, 3).**

- **The flow outlet shall not be use for driving any medical device.**

Note also that products colour (specially flow control knob) might not follow any gas colour coding.

6. Operations

6.1 Before use

Visual Inspection before fitting

- Check if there is visible external damage to the product (including product labels and marking). If it shows signs of external damages, remove it from service and identify its status.
- Visually check if the product is contaminated; if needed, for the regulator, use the cleaning procedure detailed in this section.
- Check if the service or the disposal time has not been exceeded, using GCE or owner's date coding system. If service or disposal time has been exceeded, remove the low pressure regulator from service & suitably identify its status.
- Set flow control knob on ZERO position - until the flow control knob engages.

- **Remove caps from inlet and/or flow outlet. Keep caps in a safe place for reuse during transport or storage.**

- **The low pressure regulator is dedicated only for use with the gas specified on its labelling. Never try to use for another gas.**

Fitting to medical gas pipeline system terminal units

The low pressure regulator shall be connected only at specified position. Zero flow ("0") on the flow control knob shall be positioned up (at 12'clock).

Leakage test

- After connection to medical gas pipeline system terminal unit or to medical regulator, visually check possible leakages:
 - low pressure regulator inlet connection to terminal unit or to medical regulator,
 - inlet connector to low pressure regulator body,
 - user outlet.

- **If any leakage is detected, use the procedure in chapter 6.3 and return the product to GCE for service.**

Functional test

- Check that there is a gas flow at each setting (for instance, by listening for the sound of gas flow or checking presence of bubbles in an humidifier).
- Reset flow control knob on ZERO position.

6.2 Connection & use of outlet

List of recognized accessories

To be connected to the flow outlet:

- humidifiers, breathing masks or cannulas, gas savers, nebulizers.

- **Before connecting any accessory or medical device to the low pressure regulator, always check that it is fully compatible with the product connection features & the product performances.**

Flow outlet connection

- **When connecting any accessory to the flow outlet make sure that it is not connected to the patient before operating the product.**

- Ensure the hose/humidifier is compatible with the flow outlet feature.
- Push the hose onto the regulator flow outlet/outlet of humidifier.
- Ensure the hose/humidifier is well engages.

Flow setting

- Set medical gas flow control knob on one of the available flow rates by turning the flow control knob in the requested position.

- **Always ensure that the flow control knob has engaged and not placed between two settings otherwise the flow selector will not deliver correct flow of medical gas.**

- **Do not try to apply an excessive torque on the flow control knob when it stops on the maximum flow position or in zero position.**

- **Only user and doctor is accountable for method of using at specific patient and for assessment of flow value.**

6.3 After use

- Before disconnecting low pressure regulator from the medical pipe line system terminal unit or medical regulator, reset flow control knob on ZERO position.
- Vent gas pressure from downstream equipment (i.e. via the connected consumers).
- Remove connections from user outlets.
- Refit protection caps. Before refitting the caps ensure they are clean.

7. Cleaning

Remove general contamination with a soft cloth damped in oil free oxygen compatible soap water & rinse with clean water.

Disinfections can be carried out with an alcohol-based solution (spray or wipes).

If other cleaning solutions are used, check that they are not abrasive and that they are compatible the product materials (including labels) and gas.

- **Do not use cleaning solutions containing ammonia!**
- Do not immerse in water or any liquid.
- Do not expose to high temperature (such as autoclave).

8. Maintenance

8.1 Service and Product Life time

Service

GCE recommends that a product *Periodic inspection* is undertaken every year to check proper functionality of the regulator. It should be done by skilled technician.

It should not be assumed that Periodic Inspection recommended by GCE covers every safety procedure or method required by local regulations or statutory requirements, or that abnormal or unusual circumstances may not warrant or suggest further requirements or additional procedures.

Life time

Maximum life time of the product is 10 years.

At the end of the product's life time, the product must be withdrawn from service. The owner shall put in place a relevant procedure to ensure the product cannot be used again.

8.2 Repairs

Repairs

Repairs activities cover the replacement of the following damaged or missing components:

- Quick connector
- Piston

The repairs shall be carried out by GCE authorised person only.

Any product sent back to GCE authorised person for maintenance shall be properly packaged. The purpose of the maintenance has to be clearly specified (repair). For product to be repaired a short description of fault and any reference to a claim number might be helpful.

Some repair activities concerning to the replacement of the damaged or missing components can be carried out by the owner of the product. Following parts can be replaced only:

- Caps
- Flow knob and stickers
- Hose nipple (including o-ring)

- **All labels on the equipment must be kept in good, legible condition by the owner and the user during the entire product life time.**

- **All seals and o-rings must be kept in dry, dark and clean environment by the owner and the user during the entire product life time.**

- **Use only original GCE components!**

9. Glossary

Consult operating instruction

Caution

Upper and lower humidity limit

Keep away from oil and grease

Suitable for Home care use

Suitable for Hospital care use

Suitable for Emergency care use

Keep dry

Upper and lower temperature limit

Keep away from heat and flammable material

Product serial number YY MM XXXXX
Y: manufacturing year
M: manufacturing month
X: production number

Product part number

P1

Inlet pressure range

10. Guarantee

GCE guarantees the regulator for one year or in accordance with statutory warranty rights, from date of delivery, against faulty design, material & workmanship.

GCE shall not be liable for loss of production, loss of profit or any other consequential damage or indirect loss. In the event of any fault in the goods due to defective design, materials or workmanship, our liability is limited to replacement of these goods, provided that written notification is given to GCE within three months of the date of delivery or deemed delivery, or such shorter time as may be specified in the quotation. Goods returned to GCE will not be accepted unless GCE written consent to their return has previously been obtained.

The liability of the regulator is irrevocably transferred to the owner or operator to the extent that it is modified, serviced or repaired by personnel not employed or authorized by GCE or if the apparatus is used in a manner not conforming to its intended use.

GCE cannot be held responsible for the misuse of the equipment in case of non-application of the instructions for use.

APPENDIX

Nr 1: Technical and performance data.

Nr 2: Quick coupling features and connecting / disconnecting procedure.

Manufactured by
GCE s.r.o.
Zizkova 381
583 81 Chotebor
Czech Republic

Tel : +420 569 661 111
Fax : +420 569 661 602
<http://www.gce-online.com>
© GCE s.r.o.

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