MEDIMETER®

BALL FLOWMETER

GB
CZ
HU
PL
HR

INSTRUCTION FOR USE

Document No.: 735 800 000 135
Date of issue: 01/06/2010
Revision No.: 1

1. Foreword

GCE Flowmeters are medical devices classified as class IIa according to the Medical Device Directive 93/42/EEC.

The Compliance with essential requirements of 93/42/EEC Medical Device Directive is based upon ISO 15002 standard.

2. Intended Use

Flowmeter is a device intended for the administration of the following medical gases in the treatment and care of patients:

- oxygen; - O₂
- medical air; - Air
- nitrous oxide; - N₂O

The product is not intended to be use with air or nitrogen for driving surgical tool.

This low pressure regulator is intended to be fitted to medical gas pipeline system terminal units in hospital or ambulance car or to quick connector outlet of medical regulator, with nominal pressures up to 5 bar.

3. Operational, transport and storage safety requirements

	Operations	Transport	Storage
 • Keep the product and its associated equipments away from <ul style="list-style-type: none">- heat sources (fire, cigarettes, ...),- flammable materials,- oil or grease, (especially be carefull if hand cream is used)- water,- dust.	✓	✓	✓
 • The product and its associated equipments must be prevented from falling over.	✓	✓	✓
 • Use only the product and its associated equipments in well ventilated area.	✓		
 • Always maintain oxygen cleanliness standards.	✓		

Before initial use the product should be kept in its original packaging. GCE recommends use of the original packaging (including internal sealing bag and caps) if the product is withdraw from operation (for transport, storage).

Statutory laws, rules and regulations for medical gases, accident prevention and environmental protection must be observed.

Operation conditions

 -20/+60 °C

 10/100 %

 600/1200 mbar

Storage and transport conditions

 -30/+60 °C

 10/100 %

 600/1200 mbar

- In case of storage at a temperature below -20 °C, do not operate the flow selector until it has been allowed to increase its temperature to a minimum of -20 °C.

4. Personnel instructions and training

According to Medical Devices Directive 93/42/EEC the product owner must ensure that all personnel handling the product are provided with the operating instructions & performance data and are fully trained to carry out that operation. Trainees need to be supervised by an experienced person.

5. Product description

Gas from the medical gas pipeline system is entering through the inlet connector (A). Flow is then adjusted to the requested value by turning handwheel (B), which controls the valve for flow regulation. The adjusted flow is visible on the calibrated tube via the float (C). Gas flows through the tube, when constant flow is adjusted, the float is stabilized in a certain position (a certain flow on the scale). Gas is delivered to the user through the flow outlet user connection (D).

A. Inlet quick connector


The flowmeter is fitted to the medical gas pipeline system terminal unit or to medical regulator by a gas specific male probe for specific gas. There are two basic variants: 1 – inlet connector is directly mounted on the flowmeter, 2 – connection to the rail, including hose and inlet connector.

B. Handwheel

Requested gas flow can be adjusted by the handwheel. The handwheel is connected with the valve for gas regulation. By the handwheel is also possible to close the gas flow.

C. Tube with the float

The adjusted flow is visible on the calibrated tube (l/min) via the float. When constant flow is adjusted, the float is stabilized in a certain position on the scale.

D. User outlet – outlet connection

The "flow outlet" user connection is usually with connection for hose (hose nipple) or threaded (for humidifier).

This outlet is for supplying a gas flow (l/min) atmospheric pressure directly to a patient, e.g. through a cannula or a facemask..

- **The flow outlet shall not be used for driving any medical device.**

Note also that products colour (especially handwheel) might not follow any gas colour coding.

6. Operations

6.1 Before use

Visual Inspection before fitting

- Check if there is visible external damage of the flowmeter (including product labels and marking). If it shows signs of external damages, remove it from service and identify its status.
- Visually check if the product is contaminated; if needed, for the flow selector, use the cleaning procedure detailed in this instruction for use.
- Ensure the valve for flow adjustment is completely closed (turning clockwise). When closing don't use an excessive torque (max. recommended closing torque is 0,4 Nm).
- Ensure that connected hose nipple or humidifier to the outlet has correct dimensions.

- **The flowmeter is dedicated only for use with the gas specified on its labelling. Never try to use for another gas.**

Fitting to medical gas pipeline system terminal units or to medical regulator

- The flowmeter shall be connected only in a defined position – tube up.

Functional test

- Check that there is a gas flow when turning the handwheel (for instance, by listening of the sound of gas flow or checking presence of bubbles in a humidifier).
- Close the valve for gas adjustment (turning the handwheel clockwise). When closing don't use an excessive torque.

6.2 Connection & use of outlet

List of recognized accessories

To be connected to the flow outlet:
Humidifier, breathing mask or cannula.

- **Before connecting any accessory or medical device to the flowmeter, always check that it is fully compatible with the product connection features & the product performances and there is no negative influence to outlet parameters of the flow selector.**

Flow outlet connection

- **When connecting any accessory to the flow outlet make sure that it is not connected to the patient before operating the product.**

- Ensure the hose/humidifier is compatible with the flow outlet feature.
- Push the hose onto the flowmeter flow outlet/screw onto the humidifier.
- Ensure the hose/humidifier is well engaged.

Flow setting

- Connect the flowmeter to the source of medical gas.
- Set requested medical gas flow by the handwheel (turning anticlockwise the flow increase, turning clockwise the flow decrease).

- **Do not try to apply an excessive torque on the flow control knob when it stops on the maximum flow position or in zero position.**

- **Only user and doctor is accountable for method of using at specific patient and for assessment of flow value.**

- **Fluctuating inlet pressure and ambient temperature can have influence on accuracy of flow measurement.**

6.3 After use

- Before disconnecting flowmeter from the medical pipe line system terminal unit or medical regulator, close the valve for flow regulation (turning the handwheel clockwise). When closing don't use an excessive torque (max. recommended closing torque is 0,4 Nm).
- Remove connected accessories from user outlet of the flowmeter.
- Disconnect flowmeter from the medical gas pipeline system terminal units or medical regulator.

7. Cleaning

Remove general contamination with a soft cloth damped in oil free oxygen compatible soap water & rinse with clean water. Used cleaning solution shall comply with oxygen cleanliness standards for oxygen devices.

Disinfections can be carried out with an alcohol-based solution (spray or wipes). If other cleaning solutions are used, check that they are not abrasive and that they are compatible the product materials (including labels) and gas.

- **Do not use cleaning solutions containing ammonia!**

- **Do not immerse in water or any liquid.**

- **Do not expose to high temperature (such as autoclave).**

8. Life time, maintenance and service

8.1 Product Life time

Serial number and date of production

Form of nine digit serial number stamped on the product is following:

YY MM XXXXX

YY: year of production

MM: month of production

XXXXX : sequence number

Example: serial number 050300521 shows the regulator produced in March 2005, with sequence number 521.

Product life time and waste management

Maximum life time of the product is 10 years.

At the end of the product's life time, the product must be withdrawn from service (max.10 years). The owner shall put in place a relevant procedure to ensure the product cannot be used again and follow "Directive 2006/12/EC of the European Parliament and of the Council of 5 April 2006 on waste".

8.2 Service, repairs and maintenance

Service

GCE recommends that a product periodic inspection is undertaken at least once a year to do leakage test and check proper functionality of the flowmeter.

Leakage test

- After connection to medical gas pipeline system terminal unit or to medical regulator, visually check possible leakages:

- leakage test of the valve for flow adjustment:

1. close the valve for flow adjustment (turning the handwheel clockwise). When closing don't use excessive torque (max. recommended closing torque is 0,4 Nm).
2. connect the flow outlet and hose nipple with a hose, the end of the hose immerse into a jar filled with water.
3. connect the flowmeter to the source of compressed air max. 5 bar.
4. the test: leave the flow selector 1min under pressure. Number of bubbles shall not be higher than 3 bubbles per 1 minut.
5. when test ready, disconnect the flowmeter from the source of compressed air.

- external leakage test:

1. close the valve for flow adjustment (turning the handwheel clockwise). When closing don't use excessive torque (max. recommended closing torque is 0,4 Nm).
2. close the outlet connection of the flowmeter using a plug, the plug has to have the same thread as the outlet connector of the flowmeter,
3. connect the flowmeter to the source of compressed air max. 5 bar.
4. open slowly the valve for flow adjustment (turning the handwheel anticlockwise).

5. possible leakage you will detect by application of soap water on product joints.
6. after test disconnect the flowmeter from the source of compressed air, remove the plug from outlet connection and carefully clean the surface of the flowmeter, see chapter 7 (cleaning).

- **If any leakage is detected, use the procedure in chapter 6.3 and return the product to GCE for service.**

- **The product shall not be disassembled or repaired when connected to the source of the pressure.**

Repairs activities cover the replacement of the following damaged or missing components:

- quick connector (probe),
- closing mechanism
- flow head

The repairs shall be carried out by a GCE authorised person only.

Any product sent back to a GCE authorised person for maintenance shall be properly packaged. The purpose of the maintenance has to be clearly specified (repair, maintenance). For product to be repaired a short description of fault and any reference to a claim number might be helpful.

Some repair activities concerning to the replacement of the damaged or missing components can be carried out by the owner of the product. The following parts can be replaced only:

- Caps,
- Hose nipple (including o-ring),
- Stickers,
- Handwheel

- **All labels on the equipment must be kept in good, legible condition by the owner and the user during the entire product life time.**

- **Use only original GCE components!**

9. Glossary

Consult instruction for use

Caution

Keep away from heat and flammable material

Keep away from oil and grease

Upper and lower humidity limit

Suitable for Home care use

Suitable for Hospital care use


Suitable for Emergency care use

Product serial number

Product number

	Upper and lower temperature limit		Batch number
	Keep dry		Fragile
	Use by date		Flow reading on the upper edge of the float
	Atmospheric pressure		Date of production
	Manufacturer		Outlet parameter
	Inlet parameter	P1	Inlet pressure

10. Warranty

GCE guarantees the flowmeter for one year or in accordance with statutory warranty rights, from date of delivery, against faulty design, material & workmanship. GCE shall not be liable for loss of production, loss of profit or any other consequential damage or indirect loss. In the event of any fault in the goods due to defective design, materials or workmanship, our liability is limited to replacement of these goods, provided that written notification is given to GCE within three months of the date of delivery or deemed delivery, or such shorter time as may be specified in the quotation. Goods returned to GCE will not be accepted unless GCE written consent to their return has previously been obtained.

The liability of the flowmeter is irrevocably transferred to the owner or operator to the extent that it is modified, serviced or repaired by personnel not employed or authorised by GCE or if the apparatus is used in a manner not conforming to its intended use.

GCE cannot be held responsible for the misuse of the equipment in case of non-application of the instructions for use.

APPENDIX

Nr.2: Quick coupling feature and connecting/disconnecting procedure.

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© GCE s.r.o.